

3A assembly can combine two Biobricks together into a new DNA vector.

Materials:

Linearized Plasmid Backbone (with a different resistance to the plasmid backbones containing your part samples)

Part 1 , in a different backbone than the linearized plasmid backbone

Part 2 , in a different backbone than the linearized plasmid backbone

ddH₂O

10X Digest buffer

Digest Enzyme PstI

Digest Enzyme EcoRI

Digest Enzyme DpnI

Digest Enzyme SpeI (BcuI)

Digest Enzyme XbaI

10X T4 DNA Ligase buffer

T4 DNA Ligase

BSA

Procedure:

1、 Enzyme Master Mix for Plasmid Backbone (25ul total, for 5 rxns):

5 ul NEB Buffer 2

0.5 ul BSA

0.5 ul EcoRI

0.5 ul PstI

0.5 ul DpnI (Used to digest any template DNA from production)

18 ul ddH₂O

2、 Enzyme Master Mix for Part 1 (25ul total, for 5 rxns):

5 ul NEB Buffer 2

0.5 ul BSA

0.5 ul EcoRI

0.5 ul SpeI

18.5 ul ddH₂O

3、 Enzyme Master Mix for Part 2 (25ul total, for 5 rxns):

5 ul NEB Buffer 2

0.5 ul BSA

0.5 ul XbaI

0.5 ul PstI

18.5 ul ddH₂O

4、 Digest Plasmid Backbone: Add 4 ul linearized plasmid backbone (25ng/ul for 100ng total) Add 4 ul of Enzyme Master Mix;

Digest Part 1: Add 4 ul Part A (25ng/ul for 100ng total); Add 4 ul of Enzyme Master Mix;

Digest Part 2: Add 4 ul Part B (25ng/ul for 100ng total); Add 4 ul of Enzyme Master Mix;

Digest all three reactions at 37C/30 min, heat kill 80C/20 min.

- 5、 Add 2ul of digested Plasmid Backbone (25 ng)
- 6、 Add equimolar amount of Part 1 (EcoRI SpeI digested) fragment (< 3 ul)
- 7、 Add equimolar amount of Part B (XbaI PstI digested fragment) (< 3 ul)
- 8、 Add 1 ul T4 DNA ligase buffer.
- 9、 Add 0.5 ul T4 DNA ligase
- 10、 Add water to 10 ul
- 11、 Ligate 16C/30 min, heat kill 80C/20 min
- 12、 Transform with 1-2 ul of product